

BlistAssist™



BlistAssist

An innovative medication blister card packaging assistant that brings accuracy, traceability and security to a manual process

Opportunities to capture



Increased accuracy



Improved traceability



Protection of patients

BlistAssist

Engineered to meet safety and traceability needs, it addresses all the requirements of a pharmacy's operation



A revolutionary approach to improving the manual preparation and verification of medication blister cards through a software driven solution



SureMed card prepared by the BlistAssist

BlistAssist Benefits

Accuracy

- Features a light guided tray to facilitate the tedious, time consuming and labour intensive manual preparation process
- Barcode and pill image driven process which greatly improves accuracy and reduces potential errors
- A software monitors preparation progress, allowing the user to resume production after an environmental interruption

Traceability

- The unparalleled ability to track lot numbers and expiry dates in manually prepared blister cards
- Operator traceability with a log-in and password


Safety

- Using an innovative layer verification process creates a quick and convenient verification pattern
- In addition to ensuring that no blisters are missed, it captures final preparation pictures for archiving purposes
- Provides a unique remote verification feature which liberates the preparation area and promotes workflow flexibility

BlistAssist Specifications

Approximate Dimensions	
(W)	23"
(H)	26"
(D)	22.75"
Work Stations	
All-in-one PC with 2 card set-up	
Scalability	
Additional scalable device can be connected to the primary system	

The unique exploded pill-by-pill view simplifies and accelerates verification by displaying each cup's components



Through our wide array of innovative pharmacy automation solutions and technologies, **McKesson Canada Automation** is helping improve patient care and safety.

McKesson Canada is made up of many businesses all serving the health care industry.

Our solutions enable health care service providers to get closer to the millions of patients they serve every single day, while contributing to the quality and safety of care in Canada.

McKesson Canada Automation
1.888.606.6337 extension 4
automation.sales@mckesson.ca
mckesson.ca/automation

